

Exporting results files

Introduction	This topic describes how to export results data files.
Description of results files	The software generates results files in comma-separated value (CSV) format, which is readable by a spreadsheet application such as Microsoft Excel®. You can export baseline results and performance results.
Procedure	To export results files: 1. In the Reports tab, select the results file you want to export. 2. Select File > Export. 3. In the Export dialog, enter a file name. 4. Select a destination folder, then click Save. The file is saved in CSV format to the specified folder.
More information	• About cytometer configurations (page 69) • About CS&T reports (page 54)

Setting up an experiment for offline analysis

Introduction	This topic describes how to set up experiments on a workstation that is not connected to a cytometer. The setup process is first done on an acquisition workstation and then completed on an analysis workstation. Make sure that the exported cytometer configuration name is a unique name that does not already exist on the analysis workstation. Modify the name before importing the configuration if needed.
--------------	---